

APPENDIX 8Analytical techniques used by Comalco, Andel.A. COMALCO LABORATORIES (COMALCO RESEARCH CENTRE).1. STREAM SEDIMENT AND SOIL SAMPLES.Cu, Pb, Zn, Bi, Mo, Co, Ni:

Determination by AAS after dissolution of 1 g sample in hot HCl + HNO<sub>3</sub> (Scheme 1).

Aq, Mo, Bi:

Determination by AAS after dissolution of 5 g sample in hot HNO<sub>3</sub>. (Scheme 2).

Fluorine: (total)

Determination by specific ion electrode after fusion of 0.25 g sample with flux of Na<sub>2</sub>CO<sub>3</sub>, KNO<sub>3</sub>, K<sub>2</sub>CO<sub>3</sub>, SiO<sub>2</sub> and water leach. (Scheme 8).

2. ROCK SAMPLES.Cu, Pb, Zn, Bi, Mo, Co, Ni, Cr:

Determination by AAS after attack of 1 g sample in hot HCl + HNO<sub>3</sub> + HClO<sub>4</sub>, evaporated to fuming HClO<sub>4</sub>, and dissolution in HCl. (Scheme 3). Where silicates are suspected the method is the same but attack is by HF, HNO<sub>3</sub> + HClO<sub>4</sub>. (Scheme 5).

Aq, Mo, Bi, Sb:

Determination by AAS after dissolution of 5 g sample in hot HCl + HNO<sub>3</sub>. (Scheme 4).

W:

Fusion of 0.25 g sample with KHSO<sub>4</sub>, leach with HCl, add SnCl<sub>2</sub> solution, determined with Zn-dithiol indicator in a spectrophotometer, compared with known solutions. (Scheme 7).

Au:

10 g of sample roasted, leached to dryness with HCl + HNO<sub>3</sub>, dissolved in HCl, transfer to organic solvent and determine by AAS.

Fluorine:

As for stream sediment and soil samples.