



Our reference : BU015727
 Your reference : Batch 200199
 Project code : Magnesite Project
 Report date : 11/02/99
 Report status : Final
 Page : 1 of 6

Analabs Pty. Ltd.
 ACN 004 591 664
 14 Thirkell St, Burnie
 Tasmania 7320
 Telephone : (03) 6431 6837
 Facsimile : (03) 6431 8890

ANALYTICAL DATA

Sample	MgO	CaO	SiO ₂	Fe ₂ O ₃	
MC38 126.0-127.0	43.04	1.66	0.90	4.08	120.0-130m.
MC38 127.0-128.0	42.52	2.88	1.75	2.98	
MC38 128.0-129.0	42.09	4.12	0.93	2.54	
MC38 129.0-130.0	41.98	3.70	1.47	2.50	
MC38 130.0-131.0	43.71	3.23	<0.05	2.02	
MC38 131.0-132.0	43.52	3.47	<0.05	1.90	almost
MC38 132.0-133.0	42.15	4.98	0.16	1.91	
MC38 133.0-134.0	42.45	4.28	0.40	1.94	
MC38 134.0-135.0	41.78	5.14	0.63	1.98	
MC38 135.0-136.0	42.75	2.93	2.01	1.76	
MC38 136.0-137.0	43.45	2.06	2.61	1.62	135-141 6m. 42.84MgO. 2.50CaO
MC38 137.0-138.0	43.12	2.21	2.93	1.79	
MC38 138.0-139.0	42.09	3.11	4.02	1.65	
MC38 139.0-140.0	42.66	2.37	3.53	1.76	
MC38 140.0-141.0	43.00	2.33	3.29	1.66	
MC38 141.0-142.0	42.06	3.28	3.22	1.53	
MC38 142.0-143.0	41.47	3.80	3.81	1.70	
MC38 143.0-144.0	40.23	3.87	6.04	1.77	
MC38 144.0-145.0	41.98	3.38	2.76	1.99	
MC38 145.0-146.0	41.29	3.94	2.78	2.30	
MC38 146.0-147.0	42.36	3.05	2.03	2.30	
MC38 147.0-148.5	42.81	2.76	1.12	2.98	
MC38 154.5-155.5	40.91	4.48	3.77	1.77	
MC38 155.5-156.5	40.06	4.73	4.27	2.00	
MC38 156.5-157.5	41.23	4.42	3.06	1.84	
MC38 157.5-158.5	42.07	3.97	2.42	1.57	almost
MC38 158.5-159.5	42.41	3.97	1.91	1.79	
MC38 159.5-160.5	42.15	4.22	1.50	1.70	
MC38 160.5-161.5	40.36	5.56	3.42	1.81	
MC38 161.5-162.5	41.63	4.23	2.20	1.80	
MC38 162.5-163.5	41.51	4.24	2.05	1.97	
MC38 163.5-164.5	41.53	4.02	2.78	2.06	
MC38 164.5-165.5	41.08	4.71	2.00	1.99	
MC38 165.5-166.5	41.78	3.76	2.64	2.18	
MC38 166.5-167.5	41.25	4.00	2.47	2.27	
MC38 167.5-168.5	40.58	3.50	4.31	2.54	
MC38 186.7-188.0	42.41	4.94	0.97	1.28	
MC38 188.0-189.0	42.79	3.09	5.26	1.16	
MC38 189.0-190.0	44.62	1.64	3.40	1.20	
MC38 190.0-191.0	45.27	1.56	1.20	1.17	
MC38 191.0-192.0	44.95	2.07	1.37	1.07	189-200m 11.0-44.44MgO. 2.36CaO
MC38 192.0-193.0	45.03	1.41	3.08	1.08	
MC38 193.0-194.0	44.38	1.98	3.26	1.01	
MC38 194.0-195.0	43.87	2.91	2.65	1.02	
MC38 195.0-196.0	43.96	3.32	1.55	0.98	
MC38 196.0-197.0	44.47	2.96	1.34	0.97	
MC38 197.0-198.0	43.89	2.91	2.02	0.94	
MC38 198.0-199.0	44.14	2.49	3.39	1.01	
MC38 199.0-200.0	44.31	2.70	1.60	0.97	
MC38 200.0-201.0	43.47	3.91	1.49	0.93	
Method	X408	X408	X408	X408	
Units	%	%	%	%	
Detection Limit	0.01	0.01	0.05	0.01	

Notes: N.A. = not analysed, -- = element not determined, I.S. = insufficient sample, L.N.R. = listed not received