



LEGEND:
 ● Power auger
 ○ Hand auger
 Grid position

30 ppm Zn original sample
 15 ppm Zn repeat and check duplicate sample
 30 ppm Zn duplicate sample

ANALYTICAL METHODS:
 Zn by AAS following hot conc.
 HCL and HCL/HNO₃ leach in other stages for 1 hour of 0.25g sample.

CONTOUR INTERVAL:
 19-30 50-100ppm 19-22 400-800ppm
 19-10 100-200ppm 19-50 > 800ppm
 19-15 200-400ppm

NOTE:
 This map has been compiled from enlargements of the Terrestrial Leach Department's 1:12500 scale topographic machine plots.

209253

DATE: JUNE '78
 GEOLOGIST: C.D.S.
 DRAWN: J.P.M.
 CHECKED: C.B.

GEOPEKO LIMITED 78-1313
 KING ISLAND GROUP
 No. KT27/76 V9-5
 E.L. 27/76 ELLIOTT BAY, TASMANIA 052
 VOYAGER 9
 C-HORIZON GEOCHEMICAL RESULTS
 ZINC CONTOURS

