

Test PC2/F7 - Stage grind to 100% -150  $\mu$ m at 60% solids (4 min., 500 g charge) sulphide pre-float; stage tin rougher floats.

(a) Flotation Conditions and Reagents

Stage	Conditions			Reagent Addition kg/t			
	Time (min.)		pH	PAX	PTAA	H <sub>2</sub> SO <sub>4</sub>	MIBC
	Conditioning	Flotation					
Sulphide Ro Float 1	2	5	5.7	0.01	-	-	0.04
Sulphide Ro Float 2	2	5	5.7	0.01	-	-	0.02
Tin Ro Float 1	5	5	5.0	-	0.25	0.084	0.02
Tin Ro Float 2	5	5	5.0	-	0.25	0.066	0.01
Tin Ro Float 3	5	5	5.0	-	0.20	0.054	0.01
Tin Ro Float 4	5	5	5.0	-	0.20	0.052	0.01
Tin Ro Float 5	5	5	5.0	-	0.20	0.052	0.02
Tin Ro Float 6	5	5	5.0	-	0.20	0.052	0.01
Tin Ro Float 7	5	5	5.0	-	0.20	0.048	0.01
Tin Ro Float 8	5	5	5.0	-	0.20	0.050	0.01
Tin Ro Float 9	5	5	5.0	-	0.20	0.050	0.01
Tin Ro Float 10	5	5	5.0	-	0.20	0.052	0.01
Tin Ro Float 11	5	5	5.0	-	0.20	0.052	0.01
Tin Ro Float 12	5	5	5.0	-	0.20	0.050	0.01