

Sample No.	-80# Total Extraction						-80# Cold Extraction					
	Cu	Pb	Zn	Mn	Co	Fe%	Cu	Pb	Zn	Mn	Co	Fe%
21372 152 ✓	40	40	70	150	20	2.3	20	-	30	80	-	0.2
21373 153 ✓	60	80	140	3200	50	2.2	50	60	90	1700	10	0.7
21374 154 ✓	30	30	10	40	-	0.3	20	10	-	20	-	0.2
21375 155 ✓	30	30	50	430	20	1.6	20	10	20	320	-	0.5
21376 156 ✓	100	160	90	5000	110	4.1	50	120	30	2300	30	0.9
21377 157 ✓	210	40	120	380	40	7.6	70	10	60	220	10	1.4
21378 158 ✓	80	60	90	230	30	2.9	30	40	10	120	-	0.8
21379 159 ✓	120	60	190	1500	70	5.5	30	20	50	1100	20	1.1
21380 160 ✓	200	70	200	2500	100	6.0	70	40	50	1500	30	0.8
21381 161 ✓	190	100	150	2400	90	10.0	70	20	60	1800	40	3.0
21382 162 ✓	130	180	240	3400	100	6.4	40	140	110	2100	40	1.9
21383 163 ✓	140	40	140	1500	60	8.0	40	10	30	1100	20	1.7
21384 164 ✓	70	30	70	120	20	2.6	30	10	20	50	-	0.8
21385 165 ✓	90	60	160	1000	30	5.4	40	30	70	750	-	1.0
21386 166 ✓	10	30	30	60	-	0.6	10	-	10	30	-	0.1
21387 167 ✓	10	50	50	120	-	0.6	10	20	30	90	-	0.2
21388 168 ✓	70	90	240	2800	-	1.8	50	40	150	1400	-	0.2
21389 169 ✓	50	120	100	770	20	1.8	20	70	40	420	-	0.4
21390 170 ✓	90	70	190	1800	60	5.4	30	20	70	950	10	0.7
21391 171 ✓	10	40	40	200	-	0.5	10	-	30	150	-	0.2
21392 172 ✓	10	50	40	400	-	0.8	-	10	30	270	-	0.4
21393 173 ✓	-	10	-	-	-	-	-	-	-	-	-	-
21394 174 ✓	-	20	20	20	-	0.2	-	10	10	10	-	-
21395 175 ✓	-	20	20	40	-	0.4	-	-	10	30	-	-
21396 176 ✓	-	10	-	-	-	0.1	-	-	-	-	-	-
21397 177 ✓	80	60	60	250	10	5.0	30	-	10	80	-	0.4
21398 178 ✓	60	60	70	480	20	3.8	30	10	10	320	-	0.4
21399 179 ✓	-	30	20	60	-	0.5	-	-	-	40	-	-
21400 180 ✓	40	70	160	1100	20	2.5	30	20	70	600	-	0.3
21401 181 ✓	40	50	20	70	-	1.9	30	-	-	20	-	0.4
21402 182 ✓	60	80	50	100	-	3.2	40	10	-	10	-	0.2
21403 183 ✓	-	30	30	260	-	0.5	-	-	10	180	-	0.2
21404 184 ✓	-	20	-	-	-	0.2	-	-	-	-	-	-
21405 185 ✓	-	30	10	50	-	0.2	-	-	-	30	-	-
Detection Limit	10	10	10	10	10	0.1	10	10	10	10	10	0.1